

ADATTO™

Hospital Bed

Simplicity and comfort

PATIENT WARD



New Adatto hospital bed for simplicity and comfort in patient care

Merivaara's new Adatto is a cost-effective and long-lasting hospital bed, suitable for different needs of modern patient care. It offers convenience and comfort for both caregivers and patients.

- Economic
- Ergonomic
- Easy to use



For caregivers, the Adatto hospital bed offers good working ergonomics. Its wide height adjustment range and variable treatment positions facilitates good care.

Thanks to its lightweight structure and directional castor it is also easy to manoeuvre and move around.

Patient comfort is enhanced through the 4-sectional mattress base. Additionally, back section slides when adjusted thus preventing the patient from sliding off.

Additionally, patients can adjust the height themselves with the simple hand control, making it easier for them to get into or out of the bed. The hand control adjustments can also be locked (optional feature), if needed, for increased patient safety.

Adatto complies with the new international standard for medical beds IEC 60601-2-52.

Features

- 4-sectional mattress base with integrated length extension and sliding back section for excellent ergonomics
- Ventilated cover plates made of durable and environmentally friendly, recycled ABS plastic
- Fully electric height, back section and leg section adjustments
- Solid and stable structure, designed for long-term use
- Cost-effective choice with minor investment and low maintenance costs
- Safe working load 250 kg



As a safety feature, the otherwise fully electric bed has gas-spring assisted functions for Trendelenburg adjustment and quick release of the back section.



Adatto also offers the possibility to adjust the lower leg section by a simple manual function.



Adjust the height, and the back and leg sections with an easy-to-use hand control. Adjustments can be locked, if needed (optional feature).



The 4-sectional mattress base offers variable treatment positions which are of the utmost importance for patient well-being. Back section slides when adjusted thus preventing the patient from sliding off.



Trendelenburg and reverse Trendelenburg are quick to adjust manually, when needed.

Accessories




Infusion bag holder for lifting pole
100000713



Infusion rod with
2-hooks 100005780
4-hooks 100000499



Lifting pole with triangle
A42245400



Grab handle
A41842500




Urine bottle holder
128009790



Journal basket
128009784



Monitor-/journal tray
128009768



Reading light 12V/10W - with special fixing for stretcher & bed
100000498



Battery back up
A42826400



Urine bag holder
100000799



Hand control with integrated locking function
A42826500

Technical Specifications

ADJUSTMENT RANGES

Mattress base	4-sectional with integrated length extension
Weight kg	105 kg
SWL	250 kg
Maximum patient weight	180 kg
Length (with bumpers)	2180 mm
Mattress base width	800 mm
Height adjustment	390 mm – 840 mm
Back section adjustment	0° – +71°
Leg section adjustment	-14° – + 33°
Trendelenburg adjustment	13°
Reverse Trendelenburg	4°
Castors, with central locking and directional castor	150 mm

Adjustments

Height of bed and back section	Electric
Leg section	Electric
Lower leg section	Manual
Trendelenburg and Quick Release	Manual (gas-spring assisted)

• EN 60601-1, EN 60601-1-2, IEC 60601-2-52




Solutions for efficient physical patient flow and well-being

Merivaara provides a wide range of hospital-grade furniture such as operating tables, medical lights, examination tables, trolleys and stretchers for transportation and day surgery, as well as birthing and patient beds. Additionally, our new concepts include both revolutionary integrated management systems and after-sales services.

Our versatile solutions and product offering enable convenient and cost-effective physical patient flow within the hospital, in day surgery and health clinics, as well as in nursing homes and home care. Our priority is to enhance comfort and well-being of both patients and their caregivers.

Established in 1901, Merivaara has more than 100 years of experience in designing and manufacturing hospital-grade furniture. Today our products are highly appreciated by users in more than 120 countries thanks to their ease-of-use, durability and ergonomic design.

Merivaara medical products are CE-labelled and the company's quality management system complies with EU directives for medical device. It is certified by DNV according to ISO 9001:2008 and ISO 13485:2003 standards. The company's environmental management system is ISO 14001:2004 certified.



Merivaara Corp.
Puustellintie 2, FI-15150 Lahti, FINLAND
Tel. +358 3 3394 611, fax +358 3 3394 6144
merivaara@merivaara.com
www.merivaara.com

The status of the information, specifications and illustrations in this publication is indicated by the date of publication. Whilst the greatest care has been taken over its compilation, Merivaara Corp. is not responsible for any errors or omissions. Merivaara Corp. reserves the right to make changes to the technology, features, specifications and design of the equipment and models without notice. Please check the Merivaara website for the latest updated specifications and other information.

All trademarks are Merivaara property unless otherwise stated. Patents granted or pending.